



Power Lift Recliners

AMERICA'S #1 SELLING BRAND

JUNE 2020



 GOLDEN

EZ Sleeper

Introducing the next generation of luxury power lift recliners with our patented Twilight Technology - the new EZ Sleeper with a modern, sleek design! Healthy benefits include being able to elevate the feet above the heart while experiencing a feeling of weightlessness. **Our new 5-motor system combines the ability to achieve unique positions with power headrest and power lumbar STANDARD!**



MaxiComfort EZ Sleeper

PR735 shown in Brisa® Shiitake*



Quick Ship Fabrics for PR761 & PR735

Brisa® Coffee Bean*, Brisa® Shiitake* & Brisa® Iron*



Up



TV



TV in Twilight



Lounge in Twilight



Lounge Reclined in Twilight



Zero Gravity



Rejuvenate in Twilight



Trendelenburg



Shiitake*



Coffee Bean*



Iron*

Please view fabric colors in person for accuracy.

SCAN WITH SMART PHONE CAMERA FOR "HOW-TO" VIDEO



= pocketed coil spring seat std.

*Upcharge applies

MaxiComfort SERIES



Discover Weightlessness in the MaxiComfort Cloud with our patented Twilight Technology! **Our new 5-motor system - with power headrest and power lumbar STANDARD - offers an array of rejuvenating positions to cradle the body and ease the mind.** Achieve extreme Zero-gravity and TV positions while enjoying the world's finest seating experience!

MaxiComfort Cloud with Twilight Positioning

PR515 shown in Brisa® CoffeeBean*



Quick Ship Fabrics for PR515

Hazelnut, Sterling, Brisa® Coffee Bean*, Brisa® Shiitake*

The Cloud™ power lift recliners offer our unique, ergonomic seating system for long-term comfort and support. They feature our custom bucket seat and plush backrest. Available in two styles!

SCAN FOR
MAXICOMFORT
PATENTED
POSITIONING
VIDEO



MaxiComfort Cloud
PR510-SME shown in Calypso
(Longer lead time may apply in this color.)



AutoFlex
optional.

MaxiComfort Cloud
PR510-MLA shown in Hazelnut



Quick Ship Fabrics for PR510-SME

Hazelnut, Sterling, Shiraz & Brisa® Coffee Bean*



Quick Ship Fabrics for PR510-MLA

Copper, Sterling, Hazelnut & Brisa® Coffee Bean*



= PowerPillow™ & Power Lumbar options



Hazelnut



Sterling



Shiraz



Shiitake*



Coffee Bean*



Please view fabric colors in person for accuracy.

*Upcharge applies

MaxiComfort SERIES

The DayDreamer™ is a medium/large sized power lift recliner offering our exclusive PowerPillow™ articulating headrest. It features The Cloud™ seating system and our custom bucket seat design. Recline, relax, and adjust your head and back positions with the touch of a button.



DayDreamer™
PR632 shown in Brisa® Coffee Bean*

PowerPillow™ STANDARD



SCAN FOR
PR632 VIDEO



Quick Ship for PR632M Hazelnut & Brisa® Coffee Bean*

The Relaxer, the original MaxiComfort power lift recliner is offered in two sizes, featuring four individual zipper backed pillows for truly customizable comfort.



Relaxer Medium
PR756MC shown in Sterling



Relaxer Large
PR756L shown in Copper
(Longer lead time may apply in this color.)



Quick Ship Fabrics for PR756MC & 756L: Hazelnut, Sterling & Brisa® Coffee Bean*



Hazelnut



Sterling



Coffee Bean*

Please view fabric colors in person for accuracy.



= PowerPillow™ option

*Upcharge applies

SIESTA SERIES

Siesta

Enjoy infinite positioning with our shear reducing motion seat! Power lumbar and a power headrest come standard with a USB port on the hand control. Luxurious fabrics with designer contrast stitching will stand out in any décor!



Siesta
PR445* shown in Timber

*This chair is designed by Golden and manufactured by our offshore partner.
IMPORTED PRODUCT



Power Headrest Standard



Power Lumbar Standard



Luxurious Fabric with Contemporary Contrast Stitching



Unique 2-way Motion Seat Reduces Shear



Soft-feel Hand Control with USB Charging Port



= PowerPillow™ & Power Lumbar options



Timber

Please view fabric colors in person for accuracy.



MaxiComforters

The MaxiComforters power lift recliners are now offered in FIVE SIZES, ensuring a perfect fit for just about anyone!



**America's Favorite
Lift Chair for
Over 20 Years!**

**An Old Favorite
Meets New
Upgrades in the
MaxiComforter
Power Lift Recliner!**



MaxiComforter Petite/Small
PR535PSA shown in Sandstorm



MaxiComforter Medium
PR535M shown in Brisa® Coffee Bean.
(Longer lead time may apply in this color.)



MaxiComforter Large
PR535L shown in Bittersweet.
(Longer lead time may apply in this color.)



Quick Ship for PR535-PSA:
Bittersweet, Sandstorm & Oxford



Quick Ship Fabrics for PR535M:
Bittersweet, Sandstorm & Oxford



Quick Ship for PR535L:
Bittersweet, Sandstorm & Oxford



Bittersweet



Sandstorm



Oxford



Brisa®
Coffee Bean*

Please view fabric colors in person for accuracy.



= pocketed coil spring seat std.

*Upcharge applies

SMALL TO TALL, WE BUILD THEM ALL!

Our long-standing, small-to-tall sized collection of power lift recliners now features a contemporary design and plush fabric offerings! You'll notice our popular MaxiComforter Series in a whole new light.



Refreshed Back Style Providing Maximum Comfort & Tasteful New Seam Design

Available in Luxurious Chenille Fabrics with a New Upgraded Stain Defense

Extended Arms Provide Better Utility to Enter/Exit the Chair When Lifted

New Arm Style With Additional Plush Padding for More Comfort & Value

375 - 500 lb. Weight Capacity

Golden's Most Sought-after Signature Seat, Which Provides the Utmost in Support & Comfort



MaxiComforter Large

*(shown in Anchor
color. Other colors may apply in this color.)*

MaxiComforter Tall PR535T *shown in Bittersweet*

MaxiComforter Medium Extra Wide

PR535M26 *shown in Oxford
(Longer lead time may apply in this color.)
(AutoDrive and Foot Rest Extension
not available on this model.)*


Quick Ship for PR535L:
Anchor, Sandstorm & Oxford




Quick Ship for PR535T:
*Bittersweet, Sandstorm &
Brisa® Coffee Bean**



Quick Ship for PR535M26:
Anchor, Bittersweet, Sandstorm

 = pocketed coil spring seat std.



Bittersweet




Sandstorm



Oxford



Anchor



Brisa
Coffee Bean*

Please view fabric colors in person for accuracy.



*Upcharge applies



Comforter SERIES

(a perfect fit)

The Comforter series was designed to fit individuals small to tall with weight capacities varying from 300 to 375 lbs. These chairs give the perfect sizes for the perfect fit.



Jr. Petite

PR501JP shown in Admiral

(Longer lead time may apply in this color.)



Small

PR501S shown in Pearl



Medium

PR501M shown in Palomino



Quick Ship for PR501JPT & PR501S:

Palomino & Pearl



Quick Ship Fabrics for PR501M:

Palomino & Pearl

Comforter Wides

The Comforter Wides are perfect for those who need more space between the armrests. Constructed using the highest quality components for reliable, steady lifting and reclining, these chairs offer weight capacities varying from 375 to 700 lbs.



Small 23

PR501S-23 shown in Cabernet



Medium 26

PR501M-26 shown in Pearl

(Not stocked.)



Quick Ship for PR501S-23:

Cabernet, Palomino & Admiral



Quick Ship for PR501M-26:

Palomino, Cabernet & Admiral



Palomino



Admiral



Pearl

Please view fabric colors in person for accuracy.



Large

PR501L shown in Evergreen

(Longer lead time may apply in this color.)



Tall

PR501T shown in Cabernet (Not stocked.)



Quick Ship Fabrics for PR501L:
Palomino & Pearl




Quick Ship for PR501T:
Palomino & Admiral



Did you know that Golden power lift recliners can be custom tailored to fit you perfectly? Golden Flagship Dealers are specially trained to customize a Comforter to fit you perfectly. So whether you are small, tall, or up to 700 lbs., the perfect fit is only a phone call away.



 = pocketed coil
spring seat std.

*Limited quantity. Extended lead times may apply.



SIGNATURE SERIES

(multi-functional)



**USB Port
Standard!**



Introducing the world's first three-position recliner with our patented Twilight Technology! Our dual-motor system comes standard and allows for a deeper recline to achieve rejuvenating positions no other three-position chair on the market offers.

Orion with Twilight Positioning
PR405-MLA *shown in Anchor*




Quick Ship for PR405-MLA:
Anchor, Sandstorm & Bittersweet



Anchor



Sandstorm



Bittersweet

Please view fabric colors in person for accuracy.

*Upcharge applies


INTRODUCING!



TWILIGHT Patented Seating Technology



The **ONLY**
3-position chair that
can lift your feet
above your heart!



Motor #1
Controls the movement of both the backrest and footrest.

Motor #2
Activates Twilight technology, a tilting movement that provides a weightless experience.

SCAN FOR
PR405 VIDEO



SIGNATURE SERIES

(multi-functional)

SCAN
FOR PR404
VIDEO



New compact
6-button remote



America's favorite lift and recline chair has a powered articulating headrest and power lumbar system to truly customize to your comfort from top to bottom!

Imperial PR404

shown in Brisa® Coffee Bean*

PowerPillow™ & Power Lumbar STANDARD



Quick Ship Fabrics for PR404: Hazelnut, Sterling & Brisa® Coffee Bean*

Designed for the Spine

Articulating headrest adjusts for personalized comfort

Power Lumbar



= pocketed coil spring seat std.



Hazelnut



Sterling



Coffee Bean*

Please view fabric colors in person for accuracy.



= PowerPillow™ and Power Lumbar options

*Upcharge applies



Cambridge
PR401SME
shown in Sterling



Cambridge
PR401MLA
shown in Hazelnut



Quick Ship Fabrics for PR401-SME & PR401-MLA: Hazelnut, Sterling, & Copper.



Hazelnut



Sterling



Copper

VALUE SERIES (cost-effective)

This series combines ultra-soft lumbar pillow support with our unique "grip-rite" arm design for easy transfer.



Monarch Medium
PR355M *shown in Autumn*



Capri
PR200 *shown in Cornflower*
(Longer lead time may apply in this color.)



Quick Ship Fabrics for PR200, PR355M & PR355L:
Autumn & Chestnut.



Monarch Large
PR355L *shown in Bourbon*
(Longer lead time may apply in this color.)



Autumn



Chestnut

Please view fabric colors in person for accuracy.



*Upcharge applies



Features that make Golden Power Lift Recliners America's #1 brand!

Golden Technologies is the world's leading manufacturer of power lift recliners. Founded in 1985, Golden has become America's most trusted and reliable brand. All chairs are handcrafted by American workers with domestic and imported materials (with the exception of the PR445 Siesta Series chair, which is designed by Golden and manufactured by our offshore partner). Each chair undergoes over 42 quality control inspections before the product is approved for shipping. Golden is the pioneer of innovative products and features. *Golden lift chairs are Class II Medical Devices.*



(Excluding PR445 Siesta Series. IMPORTED PRODUCT)



Self Diagnostic System



Intuitive Hand Control



Relaxing Heat and Massage



Lumbar Relief Technology



49 Pocketed Coil Seat Springs

Brisa



GOLDEN
AutoDrive™ 3.0
CUSTOMIZED COMFORT

TOTAL CONTROL

Operates the basic movement of the chair.

POSITION CONTROL

Programmable positioning buttons for individualized comfort.

CUSTOM CONTROL

Fine tunes any position by moving the backrest and footrest independently.



AutoDrive™ is the exclusive hand control system for MaxiComfort power lift recliners. This smart system has taken the guesswork out of getting into or out of your favorite positions. The AutoDrive control assists in finding favorite comfort positions and it can memorize up to four positions.

AUTODRIVE
SCAN FOR
VIDEO



* Not available on models PR535-M26 & PR502.



Cream UCR Buckskin UBS Saddle USA Bridle UBR Coffee Bean UCB Shiitake USH Iron UIR Black Onyx UBO Beet Root UBE Night Navy UNN

THE BRISA® COLLECTION



GOLDEN OFFERS THE ONLY LIFETIME* WARRANTY for all Brisa® breathable fabrics on lift chairs in the industry!

BRISA®
GOLDEN
SCAN FOR
VIDEO



** Limited Lifetime Warranty, as long as normal use is maintained. Restrictions apply. See Owner's Manual for details.*

NOTE: Crossing fabrics & series is available. Charges will apply and lead times will change.

NEW! HIGH PERFORMANCE FABRIC with ALTA™ TECHNOLOGY



For heavy-duty or high-traffic seating, Alta™ ramps up performance with a total liquid barrier, antimicrobial protection, and repellency to liquids and oil-based stains.

Exceptional cleanability comes standard with all Alta™ applications.

Whatever the space, soft and comfortable textiles are easy to maintain and tough as nails, without looking or feeling that way.



Carbon (TCA)



Limestone (TLI)



Shale (TSH)

Golden offers a 10-year warranty on High Performance Fabrics with Alta™ technology.



You Tube

www.YouTube.com/GoldenTechChannel

f Find us on Facebook

www.facebook.com/goldentechnologies



PRODUCT SPECIFICATIONS

Name	Model	Weight Capacity	Overall Width	Overall Height	Seat to Top of Back	Width Between Arms	Seat Depth	Floor to Top of Seat	Extension - While Reclined from Back of Seat to End of Footrest	Heat & Massage Option Available	Footrest Extension Option Available*	Straight Lift Option Available**	
MaxiComfort®	EZ Sleeper w/ Twilight	PR761	375 lbs.	34"	46"	30"	22"	21"	21.5"	43"	No	No	No
	Cloud w/ Twilight	PR515	375 lbs.	39.5"	46.5"	31"	20"	20"	21.5"	43"	No	No	No
	EZ Sleeper	PR735	375 lbs.	34.5"	45.5"	30"	21.5"	22.5"	20"	47"	Yes	Yes	No
	DayDreamer	PR632	375 lbs.	38.5"	45"	31"	22.5"	20.5"	20"	42.5"	Yes	No	No
	Cloud SME	PR510-SME	375 lbs.	33.5"	41"	27"	18.5"	19.5"	19"	43"	Yes	No	No
	Cloud MLA	PR510-MLA	375 lbs.	38"	45"	29"	20"	21.5"	20"	44"	Yes	No	No
	Relaxer	PR756MC	375 lbs.	36.5"	44"	29"	20"	21"	19"	43"	Yes	Yes	No
	Relaxer	PR756L	375 lbs.	39"	46"	30"	22"	22"	20.5"	44"	Yes	Yes	No
	Comforter PSA	PR535-PSA	375 lbs.	32"	40"	26"	20"	19"	20"	40"	Yes	Yes	No
	Comforter Md.	PR535M	375 lbs.	32"	42"	27"	20"	19"	20"	41"	Yes	Yes	No
	Comforter Lg.	PR535L	375 lbs.	34"	43"	29"	22"	21"	21"	42"	Yes	Yes	No
	Comforter Tall	PR535T	375 lbs.	34"	48"	31"	22"	21"	23"	45"	Yes	Yes	No
	Comforter M-26	PR535-M26	500 lbs.	37"	45"	30"	27"	20"	20"	43"	Yes	No	No
Comforter	Junior Petite	PR501JP	300 lbs.	30.5"	39"	25"	19"	18.5"	18"	38"	Yes	No	Yes
	Small	PR501S	300 lbs.	31.5"	40"	25"	20"	19.5"	19"	39"	Yes	Yes	Yes
	Medium	PR501M	375 lbs.	32.5"	42"	27"	21"	21"	19.5"	42"	Yes	Yes	Yes
	Large	PR501L	375 lbs.	34.5"	44"	28.5"	23"	22"	20.5"	43"	Yes	Yes	Yes
	Tall	PR501T	375 lbs.	34.5"	47"	30"	23"	24"	22.5"	45"	Yes	Yes	Yes
	Small X-wide	PR501S-23	375 lbs.	34.5"	40"	25"	23"	19.5"	18"	39"	Yes	Yes	Yes*
	Medium-26 Dbl.	PR501M-26D	500 lbs.	38.5"	42"	27"	27.5"	21.5"	19.5"	42"	Yes	Yes	Yes*
Signature	Orion w/ Twilight	PR405	375 lbs.	40.5"	46.5"	30.5"	23"	21"	21"	43"	No	No	No
	Imperial	PR404	375 lbs.	32.5"	43"	28.5"	19"	22"	19"	44"	No	No	Yes
	Cambridge	PR401-SME	375 lbs.	33"	41.5"	27"	20"	19"	19"	41.5"	Yes	Yes	Yes
	Cambridge	PR401-MLA	375 lbs.	35"	44.5"	31"	22"	22"	22"	44"	Yes	Yes	Yes
Siesta	Siesta	PR445	375 lbs.	34"	45"	30.5"	21.5"	19.5"	20.5"	42"	No	No	No
Value	Capri	PR200	325 lbs.	31"	41"	27"	20"	19"	18"	40"	No	No	Yes
	Monarch Md.	PR355M	375 lbs.	32"	42"	28"	21"	20.5"	18.5"	43"	Yes	No	Yes
	Monarch Lg.	PR355L	375 lbs.	34"	44"	30"	23"	21.5"	19.5"	43"	Yes	No	Yes

* May operate differently depending on model.
 ** Must be set at factory.

All chairs feature a battery backup system for power outages (9 volt batteries not included) & SmartTek® diagnostic system that communicates the electronic health of the chair.



Complies with
 UL 73 and
 CSA C22.2 No. 68





TEAM GOLDEN MAKES PENNSYLVANIA PROUD



As COVID-19 is impacting more and more individuals around the globe, Golden Technologies identified a way to use the exceptional skills of our employees in our power lift chair division and the capabilities within our manufacturing facility to quickly produce disposable face masks exclusively for a local healthcare provider. **We have refined our "IN STOCK" product offerings to meet the needs of our customers during this time.**

"I am extremely proud that we are doing our patriotic duty in defending our great nation against this virus and playing an even more important role in the communities where we live and work."

~ Rich Golden, CEO of Golden Technologies



PROUDLY MADE



BY AMERICANS

with Domestic & Imported Materials



BR-LIFEXPRESS0620

