Overall Length	33.7" / 855mm
Overall Width	16.5" / 420mm
Overall Height (at Maximum Tiller Height)	35" / 900mm
Overall Height - at Minimum Tiller Height	28" / 711mm
Legroom - Floorboard Length x Width	9.6" x 10" / 245x 252mm
Deck Height - Floorboard to Floor	6" / 150mm
Ground Clearance	1.7" / 45mm
Turning Radius	57" / 1450mm
Folded Overall Length	25" / 640mm
Folded Overall Width	16.5" / 420mm
Folded Overall Height	11" / 290mm
Weight with Std.	37 lbs / 16.7 Kg
Weight with Opt. 112V11.5Ah Li Ion Batteries	4 lbs / 1.85 Kg
Weight without Batteries	33 lbs / 14.96Kg
Seat Type	PU
Seat - Width & Depth & Height	13" x 9 " / 330mm x 230mm
Seat - Cushion Thickness	1.2" / 30mm
Seat Height - Cushion to Deck	13.4" / 340mm
Seat Height - Cushion to Floor	19.3" / 490mm
Tire Type	PU Solid
Rim Type Rear	5" / 124mm
Rim Type Front	5" / 124mm
Outside Diameter Rear	6.5"
Tire Width Rear	2" / 50mm
Outside Diameter Front	6.5"
Tire Width Front	2" / 50mm
Battery Type & Power	Lithium Ion 25.2V / 274.4WH
Brake System	Electromagnetic brake
Climbing Angle - Maximum Recommended	3 degree
Charger Type	2A off board
Charging Time	8 Hours
Controller	Dynamic 50A