



Leasure activities everywhere and for everyone !

The Hippocampe Beach & All-Terrain accompanies you through all your outdoor activities. *Versatile and comfortable,* it adapts to your activities so you can enjoy moments of *relaxation and sharing with your relatives and friends.*

Thanks to its many options for comfort and leisure, the Hippocampe becomes a true companion to reach all terrains.

Enjoy hikes in the mountains or through forests, swims and walks in the snow, *discover new horizons* !



CustOmize your Hippocampe based on your specific needs



Spoke covers

They protect the users's

hands during hiking.

Vipamat's goal is to insure a level of quality and performance that satisfies the requirements of the FDA.

ANT:



Transport bag



The transport bag protects your wheelchair during your car and plane travels.





The backrest bag allows for the transort of your personal items during your leisure activities.

Balloon wheels

user.

Protective nets

Adapt between the seat

and the armrests or the

footrest to protect the

hands and feet of the



Replace the double wheels and facilitate the wheelchair's mobility on sand for adult users.

Parking brakes

Rotating and rustproof, immobilize they the wheelchair on an incline and facilitates transfers.





Make your walks on snow easier, adapt to the terrain and prevent accumulation of snow on wheels.

Hippocampe Beach and All-terrain

Technical data



Seat length - feet S: 28.5" - 72 cm - M: 32"- 82 cm L: 36" - 92 cm - XL: 41.5" - 105 cm

Seat depth

Weight

14.5" - 37 cm



Back seat height 14.5" - 37 cm to 18.5"- 47cm


Seat height - ground

Size guide

Seat width 15.5"- 40 cm



Interior width (without armrests) 19.5"- 50 cm



Total length S : 67" - 170 cm - M : 71"- 180 cm L : 75" - 190 cm - XL : 80" - 203 cm

Weight without options : 17 kg / 37 lb

Max user weight : 130 kg / 286 lb



Total height 40.5"- 103 cm

13.5"- 34 cm

Tota 27" -

Total width 27" - 68 cm



Total length folded S : 45.5" - 116 cm - M : 49.5" - 126 cm L : 53.5"- 136 cm - XL : 58.5"- 149 cm



dele

Size S : to 4 ft 7 in - 1m40 Size M : from 4 ft 7 in - 1m40 to 5 ft 7 in - 1m70 Size L : from 5 ft 7 in - 1m70 to 6 ft 1 in - 1m85 Size XL : from 6 ft 1 in - 1m85

ᠿ

Zoom on the characteristics of the Hippocampe product range

Technical data

Seat : Comfortable hydrophobic foam (3cm thickness) covered by black fabric which dries quickly and avoids skin problems.

Frame : Stainless steel and anodized aluminum with a foam cover that is comfortale to the touch even in case of heat and cold exposure.

Push bar : Anodized aluminium tube with a handle covered in foam that improves comfort for the companion.

FDA REGISTERED

Vipamat's goal is to insure a level of quality and performance that satisfies the requirements of the FDA.

Join us !



www.vipamat.com

USA/CANADA

VIPAMAT 2316 Arbor St. Houston, TX 77004 Office : 832-516-8125 Cell : 865 385 4350 contact-usa@vipamat.com



VIPAMAT 3 Rue Gustave Eiffel 56270 PLŒMEUR - FRANCE Tél.: +33 2 97 86 24 87 Fax : +33 2 97 86 73 81 www.vipamat.com contact@vipamat.com