

Introducing our new Luggie Super, with enhanced travel scooter features including higher weight capacity and wider seat.

Compact • Foldable • Versatile



LUGGIE SUPER SPECIFICATION

LENGTH	39" (99.1cm)
WIDTH	21" (53.3cm)
WEIGHT (SCOOTER)	62 lb (28.1kg)
WEIGHT (SCOOTER+BATTERIES)	67 lb (30.4kg)
TURNING RADIUS	41.4" (104.4cm)
GROUND CLEARANCE	1.1" (2.8cm)
SEAT	Foldable Standard Seating
MAX SPEED	4mph (6.4km/h)
BATTERY	10.5 AH Lithium Battery
WEIGHT CAPACITY	360 lb (163.3kg)
WHEELS/TIRES	Front Wheels: 2" x 6" PU Solid Rear Wheels: 2" x 7" PU Solid
DRIVE SYSTEM	24-Volt DC Motor, Sealed Mini-transaxle

The information contained herein is correct at the time of publication. Factors such as battery charge, battery condition, user weight, tire pressure, and terrain type will cause variation to maximum range and speed. Due to continual product improvement and manufacturing tolerances; length, width, turning radius and ground clearance can be subject to variance of +/- 3%. We reserve the right to alter specifications without prior notice.

Luggie[®]
SUPER