

MOBIE PLUS

FULLY
ASSEMBLED

- Folds & Unfolds in Seconds
- Light Weight
- Ideal for Travel
- Fits in Any Vehicle
- Lithium Ion Battery
- Arm Rests Included
- Keyed Control Unit
- Up to 300lb Capacity
- Battery Range to 20 Miles (wth Spare Battery)
- Available in Red & Blue



MOBIE PLUS Specifications

Length:	36.9 Inches
Width:	18.1 Inches
Tires	Solid Wheels
Speed:	6 kmh/3.75 mph
Safe Gradient:	Up to 15 Degrees
Range:	Up to 15 km/9.32 miles
Turning Circle:	53.1 Inches
Weight	58.6 lbs (with battery)
Brake:	Intelligent Regenerative Electromagnetic Brake
Drive System:	Rear Wheel Drive
Handlebar:	2 Position, PU Foam
Max Capacity:	300 lbs
Seat Width:	17 Inches
Seat Height	22 Inches
Battery:	25.9V/10Ah (Lithium Battery)
Motor:	120W



Enhance Your Life

1-800-296-4151
sales@enhancemobility.com