

Elara PR118

Restful nights, energizing days.

Inspired by the moon, the DeLuna Series of luxury power lift recliners is designed to offer you the most calming comfort and restorative sleep. The Elara – the first to be introduced in this series – is medium sized with a back style that provides maximum comfort in a tasteful seam design. It features a manually adjustable three-position headrest, supportive chaise pad, and a USB port on the hand control. Available in luxurious new fabrics to match any décor. Enjoy the powerful energy and take time to better connect with yourself in your DeLuna recliner by Golden!





Chestnut Available Now!







New Fabrics Coming In Early 2021!



Crossover fabrics are unavailable.

PR118

I KI IO	Elara
Model Number	PR118
Size	Medium
Weight Capacity	375 LBS
Positioning	3
Twilight Technology	NO
Chaise Seat	YES
Pocketed Coil Seat	NO
Back Type	Seam
Removable Back Option	YES
Head Rest Cover	YES
Arm Rest Covers	NO
Overall Width	32"
Overall Height	45"
Width Between Arms	21"
Floor to Top of Seat	21"
Seat Depth	20"
Seat to Top of Back	29"
Seat to Top of Arm	7"
Distance Required from Wall Reclined	20"
Extension Reclined from Back of Seat to End of Foot Rest	42"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES
Warranty on Electrical Parts	3 Yrs.
Pro Rated Years 4 Thru 7	YES
Warranty on Mechanical Labor	3 Yrs.
Warranty on Electrical Labor	1 Yr.
Customizable*	NO
Hand Control Routing	Inside Arm
Extra Pocket Option	Standard
Heat & Massage Option	NO
Adjustable Headrest	Manual
Power Adjustable Lumbar Option	NO
AutoFlex Option Available	NO
Foot Rest Extension Option	NO
Straight Lift Option	YES
USB Charging Port	YES
Place of Origin	China

IMPORTED PRODUCT

Designed & Engineered in Old Forge, Assembled in China Using Golden Components

* See Customize Chart for Dimension limits

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



Smart Tek

In the event of a malfunction, our patented SmartTek™ diagnostic system enables you to get your lift chair up and running quickly.




Elegant seamed back design and supple urethane fabric make this chair both pleasing to the eye and highly durable.

Features & Benefits:

- Manually adjustable three-position headrest
- Detachable headrest cover
- Elegant seamed back design
- Easy grip armrest
- Soft-feel hand control
- USB port on hand control
- Extra pocket
- Contemporary plush new Imagine pattern fabric features our exclusive SofTouch™ with Stain Defense coating
- Supple urethane fabric for added softness and durability







Soft, easy grip armrests help exit the chair with ease.



Battery Backup System

Batteries not included. Battery Backup System is designed to bring the chair from the reclined to the seated position in the event of a power outage.

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.



Complies with UL 73 and CSA C22.2 No. 68

