



PRIDE[®]
M O B I L I T Y
L I V E Y O U R B E S T[®]

SIGNATURE COLLECTION



6-9

Wrangler®

Outdoor Adventures Await

10-13

Zero Turn 10

Reimagine Independence

14-17

Zero Turn 8

Navigate Tight Corners

18-21

Jazzy Air® 2

Elevate Your Life

22-25

Jazzy® Passport

Simply Fold and Travel

26-29

Go Chair®

Get Out There!



SCOTT MEUSER

CHAIRMAN & CEO
PRIDE MOBILITY®

With 33 years of experience under our belts, Pride realizes our business is about you. Our products are built with your input in mind and are designed to meet the needs of people in the real world. We offer a mobility product for every person and lifestyle.

The success of Pride Mobility Products is driven by our strong passion to make a difference in people's lives by providing innovative, market-leading products and friendly service. This begins with the dedication and loyalty of over 1,000 Pride employees: Their hard work, coupled with their skills, knowledge and expertise, supports our mission to be the friendliest company in the industry.

Together, with you, we want to continue our operational excellence by listening to your feedback and bringing it to life in our products and services. Thank you for allowing us into your homes and lives, and for being part of the Pride Family.

A handwritten signature in black ink that reads "Scott Meuser". The signature is fluid and cursive, with a long horizontal line extending to the right.





“

PRIDE SCOOTERS ARE A
GAME CHANGER!
THEY ARE SO DURABLE
AND EASY TO USE.
I HIGHLY RECOMMEND THEM!

— BARBARA L.,
SALT LAKE CITY, UT





OUTDOOR ADVENTURES AWAIT

WRANGLER®

Feel the crisp, fresh air, smell the blossoms on the trees and enjoy the peaceful serenity of the outdoors. Your adventure awaits, and nothing stands in your way. Forge a new path with the **Wrangler®**, an aggressive, outdoor scooter designed for extreme and rugged conditions. With powerful suspension, bright lighting and dual hydraulic brakes, this beast of a scooter is a powerful force.



Limited recline,
black vinyl,
high-back seat

User-friendly console

USB charger

Stone Grey

LED lights

Storage
compartment

Desert Sand

Backup sensor
with LED light

LED headlights and
directional signals

6.7 Peak HP.
dual motors

Comfort-Trac Suspension (CTS)

350 lbs.

Weight Capacity

Up To

11.4 mph

Maximum Speed

94.75"

Turning Radius

Up To

29.4

Miles Per Charge at 200 lbs.







REFINE IMAGINE INDEPENDENCE

ZERO TURN 10

Imagine your perfect scooter. A sleek, contemporary design that can traverse the most beautiful, out-of-reach places. Imagine that same scooter in your home, moving seamlessly through doors and down hallways. Your imagination becomes reality with the **Zero Turn 10**, a scooter that combines powerful dual motors and two-wheel drive with revolutionary turning technology. Let the Zero Turn 10 carry you to places unknown.

ZERO TURN 10

Candy Apple
Red

Ocean Blue

Pearl White



Our patent-pending, intelligent-turning *iTurn Technology™* ensures effortless navigation in small spaces and around corners.

400lbs.

Weight Capacity

Up To
7mph

Maximum Speed

43"

Turning Radius

Up To
24

Miles Per Charge at 200 lbs.







NAVIGATE TIGHT CORNERS

ZERO TURN 8

Say goodbye to ordinary, four-wheel scooters! It's time to turn a corner and say hello to stunning innovation and design. The Zero Turn 8 is a four-wheel scooter that offers a three-wheel turning radius. Tight corners, narrow aisles and small spaces are no match against the Zero Turn 8, which is loaded with features to make your ride comfortable and stylish.

ZERO TURN 8

- Sapphire Blue
- Pewter
- Pearl White
- Black Onyx
- Green Moonstone
- Fire Opal



Our patent-pending, intelligent-turning *iTurn Technology™* ensures effortless navigation in small spaces and around corners.

<p style="font-size: 2em; margin: 0;">325 lbs.</p> <p style="font-size: 0.8em; margin: 0;">Weight Capacity</p>	<p style="font-size: 0.8em; margin: 0;">Up To</p> <p style="font-size: 3em; margin: 0;">4 mph</p> <p style="font-size: 0.8em; margin: 0;">Maximum Speed</p>	<p style="font-size: 0.8em; margin: 0;">Up To</p> <p style="font-size: 3em; margin: 0;">38.25"</p> <p style="font-size: 0.8em; margin: 0;">Turning Radius</p>	<p style="font-size: 0.8em; margin: 0;">Up To</p> <p style="font-size: 3em; margin: 0;">13.4</p> <p style="font-size: 0.8em; margin: 0;">Miles Per Charge at 200 lbs.</p>
---	--	--	--







ELEVATE YOUR LIFE

**Jazzy
air₂**

Rise each day and see your world in a whole new way! Empower your life with the sleek and contemporary styling of the Jazzy Air® 2. With up to 12 inches of elevation in just 11 seconds, the world is within reach. Safely drive at 4 mph while elevated. No matter where your day takes you, this power chair gives you the confidence you need!

Silver
Pearl White
Ruby Red
Onyx Black
Pink Topaz
Sapphire Blue
Citrine Yellow
Garnet Red (matte)
Tanzanite (matte)
Black Pearl (matte)

Jazzy
air₂



Angle and height-adjustable headrest

Angle, height and depth-adjustable armrests

High-back or solid seat pan, reclining sport seat

Storage compartment

Drives safely at 4 mph while elevated at 12"

Integrated foot platform

Bright LED marker lights


Active-Trac® Suspension with Mid-Wheel 6® Technology

300lbs.
Weight Capacity

Up To
4mph
Maximum Speed Elevated and Non-Elevated

24.25"
Turning Radius

Up To
19.2
Miles Per Charge at 200 lbs.







SIMPLY FOLD AND TRAVEL

Jazzy[®] **PASSPORT**

It's the trip of a lifetime! Feel the anticipation as you wait for your flight in the comfort and style of your **Jazzy[®] Passport**. This compact power chair has plenty of storage and folds in just a few steps to make airport check-in effortless. When your plane touches down at your destination, you are ready to experience the adventures that await you!

Jazzy
PASSPORT



Back seat pocket

Mesh cup holder

Black

Lightweight
aluminum frame



Under seat storage

Lithium battery

Easy-access
free-wheel lever

250lbs.

Weight Capacity

Up To

3.6mph

Maximum Speed

24.25"

Turning Radius

Up To

11

Miles Per Charge at 200 lbs.







GET OUT THERE!

GOCHAIR

The night is young and nothing is stopping you. From chic restaurants to fabulous night clubs, the **Go Chair®** is your ticket to unforgettable experiences. This compact power chair is lightweight yet offers the comfort you need for hours of fun. The Go Chair® disassembles easily and offers built-in storage to keep everything you need close at hand. Embrace the nightlife and get out there with the Go Chair®!

- Ruby Red
- Pearl White
- Amber Orange
- Citrine Yellow
- Rose Quartz
- Sapphire Blue



Height and width adjustable armrests

Five-point swivel seat with latch release

Dual swing-away, under seat storage bins

One-hand feather-touch disassembly

Large footplate

300lbs.
Weight Capacity

Up To
3.7mph
Maximum Speed

Up To
25.5"
Turning Radius

Up To
13.2
Miles Per Charge at 200 lbs.



MODEL NAME / NUMBER	WRANGLER® / MV600	ZERO TURN 10 / S710ZT	ZERO TURN 8 / JAZZY ZTS
WEIGHT CAPACITY	350 lbs.	400 lbs.	325 lbs.
MAXIMUM SPEED ¹	Up to 11 mph at 350 lbs. Up to 11.4 mph at 200 lbs.	Up to 7 mph	Up to 4 mph
GROUND CLEARANCE ³	5.5" at motor	2" at motor	1.4" at motor, 1.5" at anti-tip, 3.5" (<i>side deck</i>)
TURNING RADIUS ³	94.75"	43"	38.25"
OVERALL LENGTH / WIDTH ³	61.25" / 32.5"	48" / 24.25"	45.75" / 20"
SEAT-TO-GROUND HEIGHT	30" - 31.5"	23" - 25"	21.5" - 23.5"
SEAT-TO-DECK HEIGHT	18.75" - 20.25"	18" - 20"	16" - 18"
FRONT WHEELS (<i>TIRES</i>)	14.5" tubeless pneumatic	9" solid	2" x 7" solid
REAR WHEELS (<i>TIRES</i>)	14.5" tubeless pneumatic	10.75" solid	2.5" x 9" solid
DRIVE WHEELS	N/A	N/A	N/A
SUSPENSION	Front and rear (<i>CTS Suspension</i>)	Front and rear (<i>CTS Suspension</i>)	Front and rear (<i>CTS Suspension</i>)
RANGE PER CHARGE ^{1,2}	Up to 18.6 miles (<i>standard</i>) at 350 lbs. Up to 22.6 miles (<i>standard</i>) at 200 lbs. Up to 24.2 miles (<i>optional</i>) at 350 lbs. Up to 29.45 miles (<i>optional</i>) at 200 lbs.	Up to 17.9 miles at 400 lbs. Up to 24 miles at 200 lbs.	Up to 11.6 miles at 325 lbs. Up to 13.4 miles at 200 lbs.
HEAVIEST PIECE ³	219 lbs. (<i>unit does not disassemble</i>)	75 lbs. (<i>front section</i>)	48 lbs. (<i>rear section</i>)
STANDARD SEATING	Type: Pillow-top, limited recline Weight: 48.6 lbs. Material: Black vinyl Dimensions: Width: 18.5", Depth: 16.5"	Type: Compact, foldable Weight: 30 lbs. Material: Black vinyl Dimensions: Width: 18", Depth: 17"	Type: Compact, foldable Weight: 27 lbs. Material: Black vinyl Dimensions: Width: 17", Depth: 17"
DRIVE SYSTEM	24-volt DC motor, rear-wheel drive, dual in-line	24-volt DC motor, rear-wheel drive, dual in-line	24-volt DC motor, rear-wheel drive, two-motor
DRIVE CONTROLLER	N/A	N/A	N/A
BRAKING SYSTEM	Regenerative and electromechanical	Regenerative and electromechanical	Regenerative and electromechanical
BATTERIES ⁵	(2) 75 AH (<i>standard</i>) (2) 100 AH (<i>optional</i>)	(2) 40 AH	(2) 22 AH
BATTERY WEIGHT ⁴	75 AH – 48 lbs. each; 100 AH – 60 lbs. each	30.25 lbs. each	Battery box 31.5 lbs., 13.5 lbs. each
BATTERY CHARGER	Off-board, 8 amp	Off-board, 3.5 amp	Off-board, 3.5 amp
HCPCS CODES	K0899*	K0899*	K0899*
WARRANTY	Frame: lifetime limited Drivetrain: 2-year limited Electronics: 2-year limited Battery: 6-months	Frame: lifetime limited Drivetrain: 2-year limited Electronics: 2-year limited Battery: 6-months	Frame: lifetime limited Drivetrain: 2-year limited Electronics: 2-year limited Battery: 6-months

MODEL NAME / NUMBER	JAZZY® AIR 2 / JAIR21001	JAZZY® PASSPORT / FJPT1700	GO CHAIR® / GOCHAIR
WEIGHT CAPACITY	300 lbs.	250 lbs.	300 lbs.
MAXIMUM SPEED ¹	4 mph	Up to 3.6 mph	Up to 3.7 mph
GROUND CLEARANCE ³	1"	3" at motor	1.7" at anti-tip, 1.8" at motor, 2.7" at footplate
TURNING RADIUS ³	24.25"	24.25"	25.5" with foot platform
OVERALL LENGTH / WIDTH ³	41.93" / 23.19"	40" / 23.5"	33" with foot platform / 22.5"
SEAT-TO-GROUND HEIGHT	21.25" - 33"	19"	22" - 23"
SEAT-TO-DECK HEIGHT	N/A	14.5" (<i>platform</i>)	16.25" - 17.25" (<i>platform</i>)
FRONT WHEELS (<i>TIRES</i>)	6" flat-free	8" solid casters	6" solid casters
REAR WHEELS (<i>TIRES</i>)	6" flat-free	N/A	9" solid (<i>drive</i>)
DRIVE WHEELS	10" flat-free	12" Pneumatic drive	N/A
SUSPENSION	Active-Trac® with Mid-Wheel 6® Technology	Limited	Limited
RANGE PER CHARGE ^{1,2}	Up to 16.8 miles (<i>standard</i>) Up to 19.2 miles (<i>optional</i>)	Up to 11 miles (<i>18 AH at 200 lbs.</i>) Up to 9.6 miles (<i>18 AH at 250 lbs.</i>) Up to 6.3 miles (<i>10.4 AH/travel battery at 200 lbs.</i>) Up to 5.7 miles (<i>10.4 AH/travel battery at 250 lbs.</i>)	Up to 8.7 miles (<i>standard</i>) at 300 lbs. Up to 13.2 miles at 200 lbs.
HEAVIEST PIECE ³	220 lbs. w/seat (<i>unit does not disassemble</i>)	54 lbs. (<i>unit does not disassemble</i>)	36 lbs. (<i>rear</i>)
STANDARD SEATING	High-back sport seat Material: Black vinyl Dimensions: 16" x 16" - 18"; 18" x 18" - 20"; 20" x 20" - 22"	Material: Black, Foam Dimensions: 17.5" x 16.5"	Type: Compact, foldable Material: Black vinyl Dimensions: 18" x 17" Weight: 27 lbs.
DRIVE SYSTEM	Two-motor, Mid-Wheel 6® drive	Two-motor rear-wheel drive	Two-motor rear-wheel drive
DRIVE CONTROLLER	40A, PG VR2	50A, PG VSI	40A, PG GC3
BRAKING SYSTEM	Regenerative and electromechanical	Regenerative and electromechanical	Regenerative and electromechanical
BATTERIES ⁵	(2) U1 (<i>standard</i>) (2) 40 AH (<i>optional</i>)	(1) 24V 18 AH Lithium ion (1) 25.2V 10.4 AH Lithium ion	(2) 18 AH
BATTERY WEIGHT ⁴	U1 – 24.5 lbs each; 40 AH – 32.5 lbs. each	7 lbs. (<i>18 AH</i>) 5.2 lbs. (<i>10.4 AH</i>)	14 lbs. each 29.5 lbs. battery box
BATTERY CHARGER	Off-board, 5 amp	2 amp, off-board (<i>standard, optional</i>)	2 amp, off-board
HCPCS CODES	K0899*	K0899*	K0899*
WARRANTY	5-year limited on frame Drive motors: 13-month limited Electronics: 13-month limited Battery: 13-month limited	5-year limited on frame Drive motors: 13-month limited Electronics: 13-month limited Battery: 13-month limited	5-year limited on frame Drive motors: 13-month limited Electronics: 13-month limited Battery: 13-month limited

PRIDE®

MOBILITY
LIVE YOUR BEST®

(US) 800-800-8586 | (CAN) 888-570-1113 | pridemobility.com |    



Note: Transport of the Jazzy® Passport item on commercial aircraft is subject to federal regulations and airline policy. Please consult your airline prior to your flight to ensure your Jazzy® Passport and its lithium-ion battery are permitted.

- 1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
- 2 Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 4 Battery weight may vary based on battery manufacturer.
- 5 AGM or gel cell type required.

NOTE: All specifications subject to change without notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for scooter only. Overall measurements will vary based on seating and accessory selections.

*K0899 – power mobility device, not coded by dme pdac or does not meet criteria

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

