



Smartie

The Smartie is a portable and light weight power chair designed for the user on the go. The unit has a quick-disconnect system for effortless disassembly, advanced electronics and high efficiency in-line motors. The Smartie disassembles effortlessly into three pieces, with the heaviest weighing just 29 lbs. The Smartie is also equipped with our Envirofriendly battery system, design to have a longer life span over conventional SLA batteries.

- + 250lb Weight Capacity**
- + Wireless Take-A-Part technology**
- + Adjustable Seat Height**
- + Easy to Dismantle into Three Parts**
- + Envirofriendly Battery System**



Burgundy



Blue

Model Number: UL8W

Medicare Billing Code: K0821 & E2397 (Lithium Battery)

HCPCS codes provided should not be considered as legal advice, and does not guarantee reimbursement.

Technical Specifications

Overall Dimensions (L x W x H)	in	36" x 23" x 37"
Number, Size of Tires	Front, Rear	2, 6" (2, 8")
Suggested User Weight ¹	lb	Max. 250
Weight of Heaviest Part ²	lb	29
Weight ³	lb	93
Each Battery Module Weight	lb	7.1
Battery Size	Type x PCS	Lithium Battery pack (12V12AH Equivelant)
Charger	Amp	Off board (2A)
Output Power of Motor	HP	.45 x 2
Maximum Speed	MPH	3.75
Maximum Climbing Angle/	%(°)	10(6)
Safe Climbing Angle ⁴		
Range (Per Charge) ⁶	Mile	10
Turning Radius	in	15.5"
Model #		UL8W
Standard Seat	W x D x BH	SV18171/219RG

Note: ¹ On level Road. ² Without basket and seating. ³ With basket, seat & batteries. ⁴ Base on user weight of 200LB. ⁵ Based on user weight of 165LB. ⁶ After the batteries and mechanical moving parts fully break in.