

# SmartScoot™

World's Best Travel Mobility Scooter



**SmartScoot™** is an innovative easy to use, durable, lightweight and foldable scooter for people who lead an active lifestyle, but sometimes need a little assistance getting around. Whether at airports, grocery stores, shopping malls, parks, at home or on vacation, SmartScoot™ makes your travels easier.

## **Your SmartScoot™ features include:**

### **Innovative Technology**

Fold n Lock™ Technology makes this scooter easy to fold collapse and lock in place in seconds. One of the major benefits of SmartScoot™ is that it fits easily in most automobile trunks and in all SUV's and Vans.

### **Durable Lightweight Construction**

SmartScoot™ is constructed of stainless steel for increased strength and durability, but due to its unique design weighs less than **40lbs**. To lighten your scooter for ease in lifting, simply pull off the seat and battery and the weight is reduced to **27lbs**.

### **Extended Battery Life Technology**

SmartScoot™ is powered by a removable advanced rechargeable Lithium Ion battery, which is made to travel up to **12 miles** on a single charge.

### **Made for Comfort & Stability**

With a fully adjustable steering column, seat, backrest and rear wheels, the SmartScoot™ provides maximum comfort and stability for those up to **300 lbs**.

### **Maximum Performance**

SmartScoot™ incorporates a forward and reverse mode, low, medium and high speed settings and battery strength indicator. Its weather resistant front wheel drive motor has both a disc brake and parking brake. The LED headlight and horn make travel safer and easier.

### **Airport and Airline Approved**

SmartScoot™ is approved for airline travel on most major airlines. Simply ride your SmartScoot™ through the airport and directly to the gate. When you reach your destination, your SmartScoot™ will be waiting for you at the gate.

### **Storage Solutions**

Your SmartScoot™ comes standard with a cane holder, removable luggage bar for carrying a briefcase or carry-on bag and a wire basket for holding groceries or personal belongings. A headlight, battery, and charger are also included.





# FEATURES AND SPECIFICATIONS

Specifications	Details
Scooter Size Setup	39" H x 38" L x 22" W
Scooter Size Folded with Seat	30.5" H x 38" L x 22" W
Scooter Size Folded without Seat	17" H x 38" L x 22" W
Scooter Weight Complete	39.5lbs Total 27 lbs Heaviest Part
Seat Weight	8.5 lbs
Battery Weight	4 lbs
Maximum Load Weight	300 lbs loaded or 136 kg
Scooter Speeds (3)	3, 5 and 7 mph
Ground Clearance	9" front - 3" rear
Battery	36V, 8 amp hr 288 watt Li-Ion Rechargeable
Battery Certification	UN PT. III 38.3 for Airline Travel
Battery Charger	Input Voltage : AC100V-240V
Battery Charging Time	~5 hours
Frame Materials	Stainless Steel
Brakes	Front Disk Brake
Drive	Front Wheel Drive

