

EasyTravel Mid range Scooter

Titan 4[®]



Advanced Technologies That Move You

The Titan 4W has been designed as an eminently portable and high performance scooter; in fact the perfect high performance transportable scooter. With the heaviest part weighting only 49 lbs, it is the world's most portable, adjustable & comfortable high performance scooter available. The Titan is only 24 Inches wide and has a very efficient high torque motor with the capability of climbing even steep inclines, ramps and hills.



EasyTravel Mid range Scooter

Newly designed, ergonomic, portable & foldable, indoor & outdoor electric scooter

- Allows you to keep an independent, spontaneous and active life.
- Foldable and detachable in seconds without effort or bending!
- Compact, highly maneuverable, allowing it to be transferred easily into a car, a bus, a train or a plane.
- Ergonomic location of battery pack.

Dimensions:

Length	137 cm	(54 in)
Width	62 cm	(24 in)
Height	95 cm	(37 in)
Total Weight	72 kg	(158.7 lb)
2 X 26 A/h Battery	9.5 kg	(21 lb) each
Front Frame	27 kg	(59.516 lb)
Rear Frame	26 kg	(57 lb)
Carrying Capacity	136 kg	(300 lb)
Range up to	25 km	(15.5 ml)

Weights & measurements are approximate

Available in these colors:



Kibbutz Tzora 9980300, Israel
Fax: +972-2-9991552 • Tel: +972-2-9908207
info@tzora.com • www.tzora.com

