

V6 & X8

Table of Contents

Page 2. Base
Page 2. Frame Color
Page 3. Seating System
Page 3. Reduced Drive Capability
Page 4. Power Seat Options: Seat Function Control
Page 4. Seat Sizes
Page 4. Rehab Seat Back Cane Options
Page 5. Rehab Seating Options
Page 5. Headrest Options
Page 5. Armrest Style Options
Page 6. Arm Pads
Page 6. Lateral Support Options
Page 7. Positioning Options
Page 7. Chest Harness Options
Page 7-8. Seat Restraint Options
Page 8. NonPower Legrest & Footplate Options
Page 9. Power Legrest & Footplate Options
Page 9. Tires
Page 10. Joystick Options
Page 10. Joystick Handle Type
Page 10. Joystick Mount
Page 11. Joystick Accessories
Page 11. Alternative Drive System
Page 11. Environmental Controls
Page 12. Accessories

Tech Specs, Options & Accessories

Extreme X8 & Frontier V6 Wheelchairs

Base



Front Wheel Drive



Rear Wheel Drive



X8



Mid Wheel Drive Hybrid



Mid Wheel Drive All Terrain

Color

V6 Colors



Black



Red



Mossy Green

X8 Colors



Midnight Black



Ocean Blue



Fire Engine Red



Clover Green



Tangerine



Realtree

Seating System: MUST CHOOSE ONE



2799 - No Power Seat Frame with Seating System 400 lbs. Capacity



2851 - 50° Power Tilt In Space* 340 lbs. Capacity



ASKIT - 135° Integrated Sheer Reduction (IRS) Power Recline for Rehab Seating System* 400 lbs. Capacity



1337 - 135° Power Recline Multi Positioning Seating System (MPS)* 400 lbs. Capacity



LS - Variable Seat Height with 12" Seat Elevator for Transfer Assist* 340 lbs. Capacity



2850 - Variable Seat Height & Tilt Combination* 340 lbs. Capacity

Reduced Drive Capability

RDC - This option allows for driving capability at 20% speed while elevated and/or tilt.

*Must be ordered with upgraded electronics SFCK or MFCK depending on Number of options selected

Power Seat Options: Seat Function Control



MSC5-P-S - External 5 Button Box -Placed Between the Armpad and the Joystick - controls modes through button box not through joystick



3636JS - 8 Way Toggle - Mounted on the inside or outside of joystick - Controls modes through toggle not through joystick

Seat Sizes (Width x Depth)



MPS - As small as 16" x 16" & as large as 22" x 20"



Rehab or Flat Pan- As small as 13" x 13" & as large as 23" x 21"

Rehab Back Cane Options



1367 - Standard Rehab Back Canes

Picture Coming Soon

3086 - Flip Forward Rehab Back Canes - Quick Release for Folding Back Canes Forward

Picture Coming Soon

1560 - Manual Recline Rehab Back Canes

Rehab Seating Options



CU-FN - Curve Cushion - Molded Foam w/Medial & Lateral Thigh Supports



46G - Elements Cushion - Single Layer Foam w/ Quadra Gel Pack



M2-F - M2 Cushion - Quadra Gel w/ Pelvic, Lateral & Medial Thigh Supports



BK-M - ACTA Back w/ Moldable Center Section & High Resiliency Wings - M2 Cushion-Quadra Gel w/ Pelvic, Lateral & Medial Thigh Supports



BK-R - Radius Back w/ High Resiliency Foam



DK-V - Visco Back w/ High Resiliency Foam - Back looks like an ACTA Back

Headrest Options



1370 - MPS Cloth Headrest - MPS Seat Back Only - 12x4.5x6.25



1370V - MPS Vinyl Headrest -MPS Seat Back Only - 12x4.5x6.25



HD-MF - Comfort Company Moldable Headrest Pad - 10"wide x 6" heigh *Requires HD-BH



HD-BH - Removable Multi-Axis Headrest Hardware - Depth adj .5" behind backrest to 7" forward

Armrest Style Options



1386/87 - Standard Removable Armrest - attaches to base-height adj 2"



FUA(1388/89) - Flip-Up Armrest - attaches to base-height adj 3" & angle adj 90° - 155°



PT0031 - Reclining Flip-Up Armrest - attaches to base & seat back-height adj 4" & Angle Adj 90° - 105°-



3072/3073 - Cantalever Flip Back Articulating - attaches to seat back-height 4" & angle adj 90° - 75°

Arm Pads



AP3(1377R/L) - 3"x12" Desk Length Oval



AP4(1378/79) - 4"x16" Tapered Oval



GAP LEFT(1381) - Trough Style - 6" wide x 14", 17" w/hand rest, 22" w/hand rest long



GAP RIGHT(1380) - Trough Style - 6" wide x 14", 17" w/hand rest, 22" w/hand rest long



AP75PR - Magic Mobility Leatherette Arm Pad - 3"x12" Rounded Edges



AP7550PR - Magic Mobility Leatherette Arm Pad - 3"x14" Rounded Edges



Lateral Support Options



LT1-SA - MPS Seating - Comfort Company Single Swing Away Lateral Hardware



LT1-L/RS4W4T - MPS Seating - Comfort Company Single Lateral Pad



LT1-SA - Rehab Seating - Comfort Company Single Swing Away Lateral Hardware



LT1-RGD - Rehab Seating - Comfort Company Single Rigid Lateral Hardware



LT1-RS4W4T - Rehab Seating - Comfort Company Single Lateral Pad Rigid, 5"Wx4"T, Stretch-Air Cover


Positioning Options



1509 - Removable Hip Guides - Set - Fixed Hardware, with 4" x 6" Pad



50596 - Removable Hip Guides - Set - Fixed Hardware, with 6" x 10" Pad



I2010 - Removable Knee Adductors - set, multi-axis, with 4" x 6" Pad



30094 - Flip Down Knee Adductors - set, multi-axis, flip down, with 4" x 6" Pad



Adductor Cover with Gel Insert - RQ30094

Chest Harness Options



1401 - Loop Style Positioning Chest Strap - 4" x 49.5"



SH103M - Chest Strap BODYPOINT Monoflex, Rear-Pull M - 4" x 49.5"



PBH - Padded Butterfly Positioning Chest Strap - 7" at Chest



SH331L-B1 - BODYPOINT Stayflex, Butterfly, Narrow Front-Pull L w/o Zipper - 3.5" at Chest

Seat Restraint Options



1608 - Push Button Style Seat Pos. Strap - 2"x 64"



1608XL - Push Button Style Seat Pos. Strap - 2"x 80"



ABSB - Airline Buckle Seat Pos. Strap - 2"x61"



EB205L-B1/HW312-KIT - BODYPOINT Evoflex Pelvic Stabilizer Lg & Mounting Bracker Kit - sitting hip Width 14"-18"



WC19V66 (V6)-WC 19 Transit Kit
WC19X86(X8)-WC 19 Transit Kit

Approved According to ANSI/RESNA WC 19 - Includes Lap Belt 600mm 2"x32"



WC19V67 - WC 19 Transit Kit (V6)
WC19X87 - WC 19 Transit Kit (X8)

Approved According to ANSI/RESNA WC 19 - Includes Lap Belt 700mm 2"x44"

Legrest Options

NonPower Legrest & Footplate Options



1347 - Center Mount w/1 pc fixed footplate 10" x 13", 10" height adj



1548 - Center Mount w/1 pc flip-up footplate 10" x 13", 10" height adj



1355 - Center Mount w/2 pc flip-up footplate 10" x 13", 10" height adj



1371 - Side Mount w/leg cane w/2 pc flip-up footplate 8.5" x 6", 13" height adj



AMP - Swing Away Amputee Support - 10" extension right or left



1550 - Fixed Center Mount Calf Pads - 5.5"x6.5"x1.75



SMCP - Swing Away Side Mount Calf Pads - For 1371 - w/ breathable fabric - 2.5" x 8" x 2.5"



CS - Calf Staps - Side Mount Only - Length 21"-39" Width 4.5"



SP102M - BODYPOINT Aeromech M-36" Calf Straps - Side Mount Only - Width 3.5"



SP102L - BODYPOINT Aeromech L-40" Calf Straps - Side Mount Only - Width 3.5"



1548-1R - Heel Guard - 2" height fits 1347 or 1348



827A - Foot Straps - Pair - size



HL - Heel Loops - Pair Length 11.5" -

Power Legrest & Footplate Options



1886 - Side Mount Power Elevating - Swingaway - Independent or Combined - 2"height adj



1886-CP - Side Mount Calf Pads - w/ breathable fabric - 2.5" x 8" x 2.5"



2232PG - Center Mount Power Elevating Articulating



CMCP - Fixed Center Mount Calf Pads - For 223PG - 4" x 12"



827A - Foot Straps - Pair - size



22578 - Standard 2 pc Footplate - 5"x10.5"



PT-0109 - Extra Large 2 pc Footplate - May Not Fit w/ all Seat Sizes - 5.5"x11.5"



PT-0117 - Extra Large 2 pc Footplate w/ 2" Heal Guard - May Not Fit w/ all Seat Sizes



HL - Heel Loops - Pair Length 11.5" - 13" & Width-1"

Tires



1575C - All Terrain Tire Complete w/ Brushed Aluminum Rim - Drive Tires Only



1574 - Pneumatic Grey Indoor Hybrid Tire
1573 - Solid Grey Indoor Hybrid Tire (Drive Tires Only)



1595 - Power Wheelchair Tire Jack

Joystick Options



PG-JSM - Standard Rnet Joystick Module - No LCD Display or Mono Jacks



PG-JSM509 - Color LCD Joystick Display - w/MonoJack Ports for External Switches

Joystick Handle Type



RJK - Cone-Shaped - Standard



1502 - Goal Post 4"



PC103B - Bodypoint Ball Knob - 1/4" ball



PC102A - BODYPOINT U-Shaped w/ Flex Shaft (3" or 4")



PC108A - BODYPOINT Mushroom Head Rubber Grip Handle (2" or 2.5")

Joystick Mount



1490R - Standard Fixed - Depth Adj 8.25" Right Mount

1490L - Standard Fixed - Depth Adj 8.25" Left Mount



1500P - Swing Away Joystick - Swings over 4.5" to the Right

1501P - Swing Away Joystick - Swings over 4.5" to the Left



MM - Stealth Midline Mount - 2 position mount for easy swing away - rotation 0° - 180°

Joystick Accessories



PG-AM - Attendant Control Joystick - Mounts to the back of chair or can be hand held



TBB - Tash Buddy Button - Diameter 2.5" & Height .7" Activating Force Outside edge 3.4oz - Center 7 oz



ESBLK - Egg Switch - single function push button that can scroll through modes or profiles



CPS9 - Swing Away Egg Switch Hardware - Stealth

Alternative Drive Systems



IDH300-1 - Stealth i-Drive Head Array



HDR-4M-S-C - Switch-It Head Array - Standard 3 Switch Head Array



PG-CK - Chin Control - w/ Bib Mount and Adapter



SNP - Sip and Puff - Controls the wheelchair by sipping and puffing



PG-OMNI - Omni Display - 4" Color Display for Alternative Drive Systems

Environmental Controls



PG-BTM2 - R-net Bluetooth Module 2 - Joystick controls cell phone or computer

Accessories



1605 - Waterproof Wheelchair Cover



1606 - Waterproof Joystick Cover



1603 - Medical Necessity Bag - Velcro's to Armrest 10.5" x 6" x 1.25"



1604 - Deluxe Backpack - 13.5" x 15" x 4"



PT0052 - Cup Holder - Self Leveling



A02H - Articulating Oxygen Holder



PT0035 - Cane Holder - Attaches to the seat back - May Not Work w/ Some Options

Picture Coming Soon

1396 - Removable Clear Polycarb Table - May Not Work w/ Some Options



2997(V6) / 2967(X8)- Luggage Rack 10.5" x 13"



2855(V6) / 2857(X8) - Lighting Kit Mounted on Seat - 2 Headlights, 2 Red Marker Lights w/Turn Signals



PT0029 - Set of Mud Guards - AT Drive Tires Only



PT0036(V6) / PT0037(X8) - Lift Brackets - Used with 4 Pt Boom Lift



PWB - Power Buddy - Includes USB Charging Port & Reading Light



PT0133 - Cell Phone Holder - Adj. to fit most cell phones



WCR630 - 6' Multifolding Ramp - Traction Tape



CM - Accessory Charger - 12 Volt/car outlet & USB Port with a max draw of 6 amps



LM100 - BeAdaptive Limited Mobility Shooting Mount(for gun or crossbow)



TM100 - BeAdaptive Trigger Mechanism



FM200 - Limited Mobility Fishing Mount