





## Ventura 3 Midsize Scooter


HCPCS: K0807



INSET A



INSET B



INSET C



INSET D



INSET E

- Easily disassembles for transport
- 400 lbs. Weight Capacity
- Includes interchangeable color panels in Red and Blue
- Armrests are padded and adjustable
- Floorboard with 23" of foot space
- Ergonomic Programmable throttle control
- Easy-to-adjust delta tiller
- Height-adjustable swivel seat with fold-down backrest can be adjusted forward or backward
- Flat-free, non-marking tires and anti-tip wheels
- Rear suspension
- Headlight and large, plastic carry basket
- Easy freewheel operation

### PRODUCTS

| PRODUCTS     | Description                                  |
|--------------|--|
| VENTURA318FS | VENTURA 3-WHEEL 18 FOLDING SEAT RED/BL, 1/cs |

### SPECIFICATIONS

| SPECIFICATIONS        | Description             |
|-----------------------|-------------------------|
| HCPCS                 | K0807                   |
| Weight Capacity       | 400 lbs.                |
| Top Speed             | 6 mph                   |
| Maximum Range         | 16 miles                |
| Turning Radius        | 53"                     |
| Climbing Angle        | 7.5°                    |
| Ground Clearance      | 2.5"                    |
| Dimensions            | 48" (L) x 25" (W)       |
| Floor-to-Seat Height  | 21.5" - 23.5"           |
| Seat Type             | Folding                 |
| Seat Dimensions       | 18" (W) x 17.5" (D)     |
| Motor                 | 24V x 350W x 5100 rpm   |
| Batteries (pair)      | 12V x 33AH              |
| Charger               | 4A Offboard             |
| Brakes                | Electromagnetic         |
| Front Wheels          | 10" x 3" Flat-free      |
| Rear Wheels           | 10" x 3" Flat-free      |
| Base Weight           | 88 lbs.                 |
| Battery Weight (pair) | 28 lbs.                 |
| Seat Weight           | 26 lbs.                 |
| Heaviest Piece Weight | 60 lbs. (Front Section) |

| SPECIFICATIONS   | Description |
|--|-------------|
| Total Weight   | 142 lbs.    |
| Warranty on Frame  | Lifetime    |
| Warranty on Electronic Controller/Drive Train Components | 24 Months   |
| Warranty on Batteries                                    | 12 Months   |

| ACCESSORIES | Description |
|-------------|-------------|
| SF8010      | 1/ea        |
| SF8020      | 1/ea        |
| SF8030      | 1/ea        |
| SF8050      | 1/ea        |

| PARTS        | Description |
|--------------|-------------|
| VENT-01      |             |
| VENT-04      |             |
| VENT-05      |             |
| S31001       |             |
| LRM412103-11 |             |
| S31208B      |             |
| VENT-72      |             |
| S31122       |             |
| GT8072W38    |             |
| GT8072W43    |             |
| GT8072W41    |             |
| GT8072W42    |             |
| S31012       |             |
| GT8072W32    |             |
| GT8072W44    |             |
| VENT-76      |             |
| S41160       |             |
| VENT-03      |             |
| VENT-53      |             |
| VENT-54      |             |
| GT8072W35    |             |
| GT8072W36    |             |
| GT8072W53    |             |
| S35003-PL    |             |
| VENT-55      |             |
| VENT-61      |             |
| VENT-62      |             |
| VENT-06      |             |
| VENT-02      |             |
| VENT-07      |             |
| S31147       |             |
| VENT-88      |             |
| VENT-48      |             |
| VENT-13      |             |
| VENT-45      |             |
| VENT-50      |             |
| VENT-46      |             |
| VENT-49      |             |
| VENT-58      |             |
| VENT-59      |             |
| VENT-60      |             |
| VENT-70      |             |
| VENT-71      |             |