













Rehabilitations-Systeme AG Roggenstrasse 3 · 4665 Oftringen · Switzerland Phone +41 (0)61 487 94 44 · Fax +41 (0)61 487 94 49 · info@rehasys.ch

www.vintage-mobility.ch

Swiss industrial design

TÜV medical product approval

Cintage Mobility

Gotsley

State of the art technology









THE GATSBY WELCOME BACK TO THE FUTURE.

Vintage Mobility: perfect symbiosis between design and technology. Swiss engineering combined with the approval by TÜV as a medical product guarantee an extraordinary product with a revolutionary design.

The **Gatsby** was inspired by the classic look of the legendary Ford T model, which was still being manufactured in **1926**, the founding date of our company! With the Gatsby, Vintage Mobility has launched **the first of a whole family of vintage mobility scooters**. In 2019, we will launch the next sensation. We cannot say more at this point in time, except that we truly believe that you will be as excited as we are!

Vintage Mobility is a brand of Switzerland-based **Rehabilita**tions-Systeme AG. Our products are designed and developed in Switzerland and are manufactured in China under Swiss production management. Our quality management system complies with the highest possible standards and guarantees TÜV approval as a medical product and fulfills all self-contained national standards and requirements.

- Swiss industrial design
- TÜV approval as a medical product
- Lockable front storage
 compartment
- State-of-the-art components
- Modern gel batteries
- Captains seat, fully flexibe, newest design, classic look
- Range up to 40 km
- Speed according to national regulations: 10 18 km/h
- Motor: 950W

- Charger: 8Amp
- Spare parts
- Accessories
- Max. allowable incline: 10 degrees (18%)
- Folding armrests
- Electronic S-Drive Controller
 120 Amp
- 24 VDC motor
- Dual electronic break system
- Overall Length: 1565 mm
- Overall Width: 670 mm
- Overall Height: 1300 mm