

THE MOST COMFORTABLE RECLINER IN THE WORLD!



LEGACY



Embrace the serenity of a quiet corner in your home with a cozy, luxurious recliner. Take time for yourself with the Legacy Collection from **VivaLift!**®. Put your feet up and fully stretch out on the footrest extension. Enjoy a good novel or turn on the television. No matter what activity you choose, infinite, lay-flat positioning ensures you'll always find your most comfortable and supportive position. Plus, you'll never misplace your book or the remote with the convenient side storage pocket.

Price Mobility® (US) 800-800-8586 | (CAN) 888-570-1113 pridemobility.com | prideliftchairs.com | ❸ ② ⑤



COLLECTION	LEGACY	LEGACY
MODEL	PLR-958M	PLR-958L
POSITION	Infinite Lay Flat Relaxed Reclined Upward Lift	Infinite Lay Flat Relaxed Reclined Upward Lift
SIZE	Medium	Large
WEIGHT CAPACITY	375 lbs.	375 lbs.
SEAT-TO-FLOOR	19.5″	20.5″
SEAT DEPTH	21.5″	22″
SEAT WIDTH	22″	23″
TOP OF BACK TO SEAT	25″	26"
FULLY PADDED CHAISE	Yes	Yes
FOOTREST EXTENSION	Standard	Standard
SEATING TYPE	Total Comfort	Total Comfort
CHAIR HEIGHT	42"	43.5″
DISTANCE FROM WALL	16″	16"
OVERALL WIDTH	35″	36.5″
CHAIR WEIGHT	153 lbs.	160 lbs.
WARRANTY	7-Yr. Pro-rated	7-Yr. Pro-rated





FUNCTIONS AND FEATURES:

POWER HEADREST



POWER LUMBAR



USB REMOTE



LITHIUM BATTERY BACKUP



FOOTREST EXTENSION



POWER BACKREST

POSITION GUIDE:



RELAXED



RECLINED



UPWARD LIFT



INFINITE LAY FLAT

STANDARD FABRIC: SAVILLE SUICK SHIP







WINE

GREY

BROWN

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.



= Indicates shortest available lead time

Recommended HCPCS: E0627-Seat lift mechanism, electric, any type DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.